

## Silencer siRNA cocktail kit (RNaseIII)

For the first time you use this kit:

Add 12ml of 100% ethanol (ASC) to the 2X wash solution and check the label. From now, it will be called 1X Wash solution and store it at RT.

- (1) Bring an ice bucket.
- (2) Get the T7 Enzyme Mix out of the freezer and place it quickly on the ice.
- (3) Vortex 10X T7 Rxn buffer as well as the four NTP solutions.
- (4) Thaw the 4 NTP solutions on the ice while keep the 10X T7 Rxn buffer at RT.
- (5) Spin these microfuges prior to open them to: (a) avoid contamination and (b) prevent loss.

---

### In vitro transcription

The *in vitro* Rxn should be done at RT. This is to avoid Spremidine (present in transcription buffer) to precipitation of DNA at low temperature.

**The order of adding materials should be:**

Water + DNA template, then add 10X Rxn buffer.

**For 20 $\mu$ l transcription Rxn:**

Water up to 20 $\mu$ l

0.2-2 $\mu$ g DNA (with T7 promoter at the ends).

2 $\mu$ l of 10X T7 Rxn buffer

2 $\mu$ l of each NTP

2 $\mu$ l of T7 polymerase Mix

- (1) Flick the tube or gently pipet up and down → spin it down → incubate for 2hrs at 37C.
- (2) Reserve 0.5 $\mu$ l of this transcription Rxn for 1% agarose electrophoresis. Run 1/400 in non-denaturing agarose and dilute the sample in TE
- (3) Incubate the sample at 75C for 5minutes and cool at RT. This maximizes the duplex RNA formation. In our case where a small DNA (i.e. <800bp) and both ends have T7 promoter sequence, no need for this step. However, it is still better to do it since there will be no drawback.

**Purification of dsRNA** (by nucleic acid digestion and filtration):

Nucleic acid digestion of DNase I and RNase A

For 20 $\mu$ l of in vitro transcription tube:

-20 $\mu$ l of dsRNA

-21 $\mu$ l Nuclease free water

-5 $\mu$ l 10X digestion buffer.

-2 $\mu$ l DNaseI

-2 $\mu$ l RNase A

---

50 $\mu$ l total

Incubate not more than 2 hours. In the first 15 minute incubation, ssRNA will be digested by RNase A while any DNA can be digested in 1 hr.

**Purification step:**

To remove proteins (the enzymes), free nucleotides and oligonucleotides but retains dsRNA until it is eluted:

- (1) Preheat the elution solution to 95C
- (2) Assemble the following:
  - 50µl of dsRNA (obtained from digestion step)
  - 50µl 10X Binding buffer
  - 150µl Nuclease free water
  - 250µl 100% ethanol
- (3) Gently mix
- (4) Pipet this 500µl onto a Transcription Rxn filter cartridge
- (5) Centrifuge at max speed for 2min → Discard the flow through and place the cartridge in a new collection tube.
- (6) Add 500µl of wash solution → Centrifuge → Discard the flow through.
- (7) Repeat step (6) → Centrifuge → Discard the flow through
- (8) Centrifuge again without any wash solution.

**To recover dsRNA:**

- (1) Transfer the rxn filter cartridge to a new collection tube.
- (2) Add 50µl of **hot elution buffer** (preheated to 95C). If you use a cold elution buffer, let it stand for 2 minutes.
- (3) Centrifuge for 2 minute at highest speed.
- (4) Repeat steps 2 and 3 on the same collection tube.

**To quantitate dsRNA:**

Measure dsRNA at 260nm.

**Storage:**

dsRNA is stable at -20C when it is stored in the elution buffer.

- Run 1/400 of clean dsRNA on 2% Agarose (non-denaturing) to check up the integrity and efficiency of dsRNA formation.
- To have 1/400 of dsRNA from the original 100µl in elution buffer, do the following:
  - (1) take 1µl of 100µl dsRNA + 9µl (10mM Tris and 1mM EDTA), i.e. 1:10
  - (2) Take 2.5µl of this 1:10 + 7.5µl of TE → 1:400 dilution.

**dsRNA → siRNA (by RNase III digestion) and the purification of siRNA**

One rxn can digest up to 15µg of dsRNA/rxn. This is enough to get siRNA for about 120 transfections in a 24well dish at 50nM final concentration.

- dsRNA up to 15µg (up to 30µl)
- 15µl RNase III.
- 5µl of 10X RNase III buffer.

-to 50 $\mu$ l water (nuclease free).

(1) Incubate in water bath at 37C for 1hr

(2) Pre-wet siRNA purification unit with 50 $\mu$ l of nuclease free water, centrifuge for 8 minutes for 14,000 x g.